

## Immunofluorescence staining protocol

### Material:

- 3 µm FFPE-cuts on Advanced Adhesive StarFrost® slides
- Xylol
- Decreasing ethanol series: 100 %, 100 %, 96 %, 80 %, 70 %, 40 %, ddH<sub>2</sub>O (deionized water)
- PBS (preparation see appendix, store at room temperature)
- Plastic box (e.g. from Ikea) + lid with holes for antigen retrieval
- Citrate buffer (pH 6, mix of solution A & B, for preparation see appendix)
  - solution A: 14.4 ml
  - solution B: 65.6 ml
  - ddH<sub>2</sub>O: 720 ml
- PAP Pen
- Ice
- Humid chamber
- 10x Rotiblock (Carl Roth, A151)
- TritonX (e.g. SERVA Electrophoresis GmbH, 37240)
- Blocking/permeabilization buffer: 1x Rotiblock in ddH<sub>2</sub>O / 0.1% TritonX (fill into drop bottle, store at 4 °C)
- Staining buffer: 1x Rotiblock in ddH<sub>2</sub>O (store at 4 °C)
- 1<sup>st</sup> antibody: antibody-specific concentration
- 2<sup>nd</sup> antibody: antibody-specific concentration (for AlexaFluor fluorophor from Invitrogen 1:500)
- DAPI (Thermo Fisher Scientific, D3571)
- Fluorescence Mounting Medium (Agilent Technologies, S302380)
- Cover slips

### Method:

#### **I) Deparaffination & dehydration**

- Xylol 3 x 5-8 min
- Decreasing ethanol series: 100 %, 100 %, 96 %, 80 %, 70 %, 40 %, 1 min each
- Wash in ddH<sub>2</sub>O

#### **II) Antigen unmasking:**

- Put slide holder into plastic box with citrate buffer (always two holders à 10 glass slides per plastic box)
  - Microwave (microwave-specific settings, boiling time should be ~20 min)
  - Cool whole box on ice for at least 15 min (possible pause point)
  - Wash slides 1 x 5 min in PBS
- Encircle cuts with PAP Pen; apply blocking buffer drops immediately (do not let cuts dry out!)

Incubation with **80 - 100 µl/ cut** in humid chamber

All steps at room temperature

### III) Staining

step	reagent	time
Blocking	Blocking/permeabilization buffer	1 h
1st Antibody	Antibody-specific concentration in staining buffer	1 h at room temperature OR over night at 4 °C
Wash	PBS	2 x 5 min
2nd Antibody	In staining buffer	1 h
Wash	PBS	3 x 5 min
Counter-staining	DAPI (1:1000 in PBS)	10 min
Wash	PBS	3 x 5 min
Coverslip	Fluorescence mounting medium	
Dry	In the dark with open lid!	Over night

Afterwards store at 4 °C for up to 1 month or, for several months, at -20 °C (tape seal box air-tight)

### Appendix:

- **PBS:**

NaCl 36 g

Na<sub>2</sub>HPO<sub>4</sub> + 2 H<sub>2</sub>O 9.4 g

KH<sub>2</sub>PO<sub>4</sub> 2.15 g

add 5 l ddH<sub>2</sub>O, shake and let settle for at least 1 h

shake again before first use

- **Citrate buffer** (prepare solutions in advance and store at 4 °C)

solution A: 14.4 ml

Citric acid 21.01 g

ddH<sub>2</sub>O 1000 ml

solution B: 65.5 ml

Sodium citrate 29, 41g

ddH<sub>2</sub>O 1000 ml

ddH<sub>2</sub>O: 720 ml